



Annexin V-EGFP Apoptosis Detection Kit

Manufactured by BioVision.

ALX-850-253-KI01: ~25 tests

ALX-850-253-KI02: ~100 tests

ALX-850-253-KI03: ~400 tests

(Version 01: April 1, 2005)

NORTH AMERICA

AXXORA, LLC
6181 Cornerstone Court East
Suite 103
San Diego, CA 92121-4727
Phone: (858) 550-8828
Fax: (858) 550-8825
E-mail: axxora-usa@axxora.com

SWITZERLAND/REST OF THE WORLD

ALEXIS CORPORATION
Industriestrasse 17, Postfach
CH-4415 Lausen / Switzerland
Phone: +41 61 926 89 89
Fax: +41 61 926 89 79
E-mail: alexis-ch@alexis-corp.com

GERMANY

AXXORA DEUTSCHLAND GmbH
Marie-Curie-Strasse 8
DE-79539 Lörrach
Phone: (07621) 5500 522
Fax: (07621) 5500 523
E-mail: axxora-de@axxora.com

UK & IRELAND

AXXORA (UK) LTD.
P.O. Box 6757
Bingham, Nottingham NG13 8LS
Phone: +44 1949 836111
Fax: +44 1949 836222
E-mail: axxora-uk@axxora.com

For laboratory use only. Not for human or diagnostic use.

I. Introduction:

The **Annexin V-EGFP Apoptosis Detection Kit** is based on the observation that soon after initiating apoptosis, most cell types translocate the membrane phospholipid phosphatidylserine (PS) from the inner face of the plasma membrane to the cell surface. Once on the cell surface, PS can be easily detected by staining with an enhanced green fluorescent protein (EGFP) fusion of annexin V, a protein that has a strong natural affinity for PS. The one-step staining procedure takes only 10 minutes. In addition, the assay can be directly performed on live cells. Detection can be analyzed by flow cytometry or by fluorescence microscopy with a FITC filter. EGFP is brighter and more photo-stable than other fluorescent reagents.

II. Kit Contents:

Components	850-252-KI01	850-252-KI01	850-252-KI01
	25 tests	100 tests	400 tests
Annexin V-EGFP	125 µl	500 µl	2 ml
1x Binding Buffer	12.5 ml	50 ml	2 x 100 ml
Propidium Iodide (PI)	125 µl	500 µl	2 ml

III. Annexin V-EGFP Assay Protocol:

A. Incubation of cells with Annexin V-EGFP

1. Induce apoptosis by desired method.
2. Collect $1-5 \times 10^5$ cells by centrifugation.
3. Resuspend cells in 500 µl of 1x Binding Buffer.
4. Add 5 µl of Annexin V-EGFP and 5 µl of propidium iodide (PI, optional.)
5. Incubate at room temperature for 5 min in the dark.

Proceed to B or C below depending on method of analysis.

B. Quantification by Fluorescence Microscopy

Analyze Annexin V-EGFP binding by flow cytometry (Ex = 488 nm; Em = 530 nm) using FITC signal detector (usually FL1) and PI staining by the phycoerythrin emission signal detector (usually FL2).

For analyzing adherent cells, gently trypsinize and wash cells once with serum-containing media before incubation with Annexin V-EGFP (A.3-5).

C. Detection by Fluorescence Microscopy


1. Place the cell suspension from Step A.5 on a glass slide. Cover the cells with a glass coverslip.

For analyzing adherent cells, grow cells directly on a coverslip. Following incubation (A.5), invert coverslip on glass slide and visualize cells. The cells can also be washed and fixed in 2% formaldehyde before visualization.

Note: Cells must be incubated with Annexin V-EGFP before fixation since any cell membrane disruption can cause nonspecific binding of Annexin V to PS on the inner surface of the cell membrane.

2. Observe the cells under a fluorescence microscope using a dual filter set for FITC & rhodamine.

Cells which have bound Annexin V-EGFP will show green staining in the plasma membrane. Cells which have lost membrane integrity will show red staining (PI) throughout the nucleus and a halo of green staining (EGFP) on the cell surface (plasma membrane).



IV. Storage:

Store kit at +4°C. PROTECT FROM LIGHT!